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Cable ISE to titrator 51089954
Cable ISE to ion meter 52300004

Literature

Guideline to Ion selective measurement 51300075

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
Subject to technical changes.
02/10 © Mettler-Toledo AG. Printed in Switzerland. 52 003 845

Internet: www.mt.com

METTLER TOLEDO

English

General information

The perfectION™ comb Na⁺ sodium electrode with built-in reference is suitable for measurements of sodium activities and sodium concentrations in aqueous systems.

Specifications

Type of ion selective membrane:

glass membrane

Type of junction:

ceramic

Resistance of the membrane:

<100 MΩ

Measurement range:

1 x 10⁻⁷ mol/L to 1 mol/L sodium

Measurement temperature range:

0–80 °C

Optimum pH range:

8...11

Reference electrolyte:

3 mol/L KCl (alternatively 0.1 mol/L NH₄Cl or as described in the application brochures)

Interferences:

Silver ions (Ag⁺) must be absent. The sodium electrode is sensitive to the following ions. The ratio of this ion to sodium (X⁺/Na⁺) should thus not be larger than the value in brackets: H⁺ (<0.001), Li⁺ (<1), K⁺ (<5), NH₄⁺ (<50), Mg²⁺ (<2000).

Precision

In the absence of interfering substances, a precision of better than ±1 mV corresponding to ±4% of the measured Na⁺ concentration can be achieved.

Response times

When changing from lower to higher concentrations below 10⁻⁴ mol/L, responding takes approx. 60 seconds, above 10⁻⁴ mol/L less than 60 seconds; when changing from higher to lower concentrations several minutes.

Materials

Electrode body and membrane made of glass.

Recommended ISA solution

ISA (Ionic strength adjuster) solution for preconditioning both calibration and sample solutions. Prepare as follows: Weigh out 200 g NH₄Cl and dissolve in 500 mL deionized water, add 50 mL concentrated NH₃ solution and make up to 1 L. The ration of standard or sample to ISA should be 100:2.

Maintenance and service

Deposits on or contamination of the surface of the membrane can be removed by rinsing the sodium electrode with deionized water or a lab cleaning agent. After every treatment, the sodium electrode must be rinsed with deionized water and then conditioned for several hours in saturated sodium solution.

Storage

We recommend to store the electrode in 1 mol/L NaCl. It is an advantage to store the electrode in the protective tube supplied (protection against damage).

Conditioning

The sodium electrode needs to be conditioned in saturated NaCl solution for 48 hours.

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