

DAIHAN-brand® UV Transilluminators, “WUV”, Standard- & Compact Case-type, 48W/90W with Hinged UV-Blocking Safety Cover, Economical, Robust, Versatile, 312-/365-/312+365-nm, with Certi. & Traceability
 자외선 발산기, 기본형 핵산(DNA, RNA), 스탠다드-/컴팩트-케이스 타입



Wisdom

Transilluminator, Standard Case-type “WUV-L50”

- with Certificate & Traceability :
Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Low Cost, Robust, Versatile, Affordable & Gentle Excitation
- Hinged UV-Blocking Safety Cover Protects User from UV Radiation
- Stainless steel Frame and Robust Steel Housing
- Increased UV Intensity and Homogeneity
- CE Certified ■ PL(Product Liability) Insurance

* Ideal for Standard Gel Documentation and Preparative Work, especially for Detection of very small amount of DNA & RNA Stained by EtBr
 * Long-Life Electronic (Low Heat, no Flickering)
 * High (100%) / Low (75%) Intensity Switch

Specification

| Wave Length & Model | Single 365 nm “WUV-L10” | Single 312 nm “WUV-L20” | Dual 365+312 nm “WUV-L50” | Single 365 nm “WUV-M10” | Single 312 nm “WUV-M20” |
|---------------------------------------|--|-------------------------|---------------------------|-------------------------|-------------------------|
| Features | Standard Case-type | | | Compact Case-type | |
| UV Lamps | 15W×6ea (90W) | | | 8W×6ea (48W) | |
| Power Consumption | 95 W | | | 52 W | |
| Filter Size (mm) | 200×200 mm | | | | |
| Dimension | 495×330×h95 mm | | | 350×275×h95 mm | |
| Weight | 10 kg | | | 8 kg | |
| Packing Size (mm) & Gross Weight (kg) | 520×380×h220 mm, 12 kg | | | 380×320×h220 mm, 10 kg | |
| Power Supply* | 1 Phase, AC 120V, 60Hz or AC 230V, 50/60Hz | | | | |

* Other Specifications are available upon Customer's Request.



Wisdom

Transilluminator, Compact Case-type “WUV-M20”

Ordering Information

Standard Case-type, 15W, UV Lamp × 6ea

| | Ea. Price, ₩ |
|--|--------------|
| DH.WUV00010 Transilluminator, 230V, 365nm, “WUV-L10” | 729,900 |
| DH.WUV00020 , 312nm, “WUV-L20” | 997,500 |
| DH.WUV00050 , 365+312nm, “WUV-L50” | 871,000 |
| DH.WUV02010 Transilluminator, 120V, 365nm, “WUV-L10” | 729,900 |
| DH.WUV02020 , 312nm, “WUV-L20” | 997,500 |
| DH.WUV02050 , 365+312nm, “WUV-L50” | 871,000 |

Compact Case-type, 8W, UV Lamp × 6ea

| | Ea. Price, ₩ |
|--|--------------|
| DH.WUV00110 Transilluminator, 230V, 365nm, “WUV-M10” | 660,930 |
| DH.WUV00120 , 312nm, “WUV-M20” | 997,930 |
| DH.WUV00210 Transilluminator, 120V, 365nm, “WUV-M10” | 660,930 |
| DH.WUV00220 , 312nm, “WUV-M20” | 997,930 |

Spare UV Lamps

| | Ea. Price, ₩ |
|--|--------------|
| DH.WUV01010 UV Lamp, 365nm, 15W, for “WUV-L10/L50”, “UVL365” | 24,540 |
| DH.WUV01020 , 312nm, 15W, for “WUV-L20/L50”, “UVL312” | 79,140 |
| DH.WUV01110 , 365nm, 8W, for “WUV-M10”, “UVM365” | 22,140 |
| DH.WUV02110 , 312nm, 8W, for “WUV-M20”, “UVM312” | 73,140 |