

DAIHAN-brand® Programmable Ball Mill, “BML”, 2-/6- Pots, with Certi. & Traceability with Digital Feedback Control System, Digital LCD, 50~600 rpm, ±5 rpm, 10 Steps 6 Programs



2-Places Programmable Ball Mill “BML-2”



6-Places Programmable Ball Mill “BML-6”



Minimum Noise and Safety

Optional Soundproof Cabinet “BMC-6” for 6-Places Programmable Ball Mill “BML-6”

- with Certificate & Traceability : Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
 - Ideal for Grinding/Mixing : Crystalline Chemicals, Inks, Paints, Pigments, Plants, Racks and Soils
 - Microprocessor, Programmable, Variable Speed : 50~600 rpm
 - Digital Feedback Control and LCD : Speed (RPM), Power Bar-Graph and Timer Status
 - CE Certified
 - PL(Product Liability) Insurance
- * Antistatic Rubber Roller : Keeps Jars / Pots Centered
 - * Digital LCD with Back-Light Function
 - * Over Load Protection
 - * Ergonomic and Compact Design
 - * NBR Roller with Pot Guide Ring
 - * Programmable Operation : 10 Steps 6 Programs
 - * Storage Function (RPM, Timer)
 - * Soundproof Cabinet for Safety and Minimum Noise (Optional)

Specification

Model	2-places“BML-2”	6-places“BML-6”
Capacity	up to 2 pots	up to 6 pots
Motor	DC 40W Geared Motor (5 : 1 ratio)	DC 100W Geared Motor (5 : 1 ratio)
Maximum Load	30 kg	60 kg
Speed Range	50 ~ 600 rpm	
Speed Accuracy	±5 rpm	
Timer	99 hr 59 min (Programmable)	
Display	Digital LCD with Back -Light Function	
Controller	Programmable Digital Feedback Controller with Jog-Shuttle Switch (Turn+Push)	
Program Function	10 Steps, 6 Programs	
Safety Device	Over Load Protector for Motor	
Dimension (w×d×h)	Roller (mm)	50×300
	External (mm)	620×300×515
Packing Size (w×d×h) (mm) & Gross Weight		70×800
		1,200×500×795
	760×460×700	1,300×620×1,010
	52 kg	115 kg
Power Consumption	40 W	100 W
Power Supply & Cord/Plug*	1 Phase, AC 120 V, 60 Hz or AC 230 V, 50 / 60 Hz	

* Other Specifications are available upon Customer's Request.

Ordering Information

For 2 Pots

- DH.WBM05002 Ball Mill, 2-places“BML-2”, 230V
- DH.WBM05102 , 120V

For 6 Pots

- DH.WBM05006 Ball Mill, 6-places“BML-6”, 230V
- DH.WBM05106 , 120V

Optional Soundproof Cabinets

* Designed for Minimizing Noise, and Safety Work

- DH.WBM505102 Soundproof Cabinet, “BMC2”, for 2-places“BML-2”, w1,020×d570×h1,120mm
- DH.WBM505106 Soundproof Cabinet, “BMC6”, for 6-places“BML-6”, w1,600×d770×h1,400mm