

Cat. No	Description	Articles
---------	-------------	----------

DAIHAN-brand® Multifunction Vortex Mixer Set, "VM-10", Up to 3,300rpm, with Certi. & Traceability with ϕ 76mm(3 inch) Platform Head & Pop-off Cup Head, Variable Speed Control, Continuous or Touch-on Mixing

- with Certificate & Traceability : Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Ideal for Mixing Solutions in Test Tubes or Small Flasks with Circular, Oscillating Movement
- Variable Speed Control Allows Slow Speed Shaking Action up to High Speed Vortexing
- CE Certified
- PL(Product Liability) Insurance

- * Gentle Mixing to re-suspending Pellet
- * Ideal for Mixing Solutions in Test Tubes or Small Flasks with Circular, Oscillating Movement
- * Supplied with ϕ 76mm Platform and Pop-Off Cup
- * Realize High Power, Stable Motion with Low Noise
- * Tough Metal Housing Provides a Stable Platform for all types of Mixing
- * 3-position Power Switch Operates with Continuous Action or "Touch On" Control
- * Variable Speed Control Allows Slow Speed Shaking Action up to High Speed Vortex
- * Aluminum Die-casting with Powder Coating
- * Various Shapes, Sizes and Materials of the Head allows for mixing of almost all Common Tubes or Containers

Specification

Model	"VM-10"
Speed Range	0 ~ 3,300 rpm
Shaking Motion	Orbital
Control	Analog Phase Controller
Head	included Pop-off Cup Head and ϕ 76mm Platform
Material	Aluminum Die-casting with Powder Coating
Dimensions (w × d × h) (mm)	160 × 128 × h145 mm
Weight	3 kg
Others	CE Certified 3-position Power Switch Operates Aluminum Die-casting with Powder Coating included Pop-off Cup Head and ϕ 76mm Platform
Packing Dimension (w × d × h) (mm)	260 × 200 × h220 mm
Gross Weight	4 kg
Power Consumption	30 W
Power Supply & Cord Plug*	1 Phase, AC 120V, 60Hz or AC 230V, 50 / 60Hz

* Other Specifications are available upon Customer's Request.

Ordering Information

※ Included Acc. : "PM110" ϕ 76mm(3inch) Platform Head and "PM210" Pop-off Cup Head, Remove Applicator

- DH.WVM00010 Vortex Mixer Set, 230V, 60 Hz, "VM-10", 3,300 rpm
 DH.WVM00015 , 50 Hz, "VM-10", 3,300 rpm
 DH.WVM00020 Vortex Mixer Set, 120V, 60 Hz, "VM-10", 3,300 rpm
 DH.WVM00025 , 50 Hz, "VM-10", 3,300 rpm

Spare Heads


- DH.WVM00110 (1) Platform Head, "PM110", ϕ 76mm(3inch)
 DH.WVM00210 (2) Pop-off Cup Head, "PM210"

Optional Heads


- DH.WVM00310 (3) Microtube Platform Head(0.5~2.0ml), "PM310"



Multifunction Vortex Mixer Set "VM-10" with ϕ 76mm(3") Platform (PM110)



(1) Platform Head, ϕ 76mm(3inch) "PM110" (Included)



(2) Pop-off Cup Head "PM210" (Included)



Remove Applicator (Included)



(3) Microtube Platform Head, 0.5~2.0ml Tube "PM310" (Optional)