
Description


DAIHAN-brand[®] Premium Magnetic Stirrer, "MS-A", 180×180mm Ceramic-Coated Plate, with Certi. & Traceability with PWM Control System, Precise Speed Control, Ceramic Coated Al-Plate, up to 1,500 rpm, Max. 20 Lit.



Premium Magnetic Stirrer with Stir Bar "MS-20A"



Specialized Features



Realize Stable &

Easy to Use!

DAIHAN

196

 with Certificate & Traceability : Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
Control Plate Chamical (Acid Paristers)

- Ceramic Coated Plate-Chemical/Acid Resistance
- CE Certified
- PL(Product Liability) Insurance
- * Ideal for Heat-sensitive Procedures such as Microbiology and Biochemistry Applications
- * Mixing for Chemical, Laboratory, Pharmaceutical Materials
- * Superior Linearity for the Control of Stirring Speed by Specially Designed Control Circuitry (PWM Control)
- $* \ {\rm Smart} \ {\rm Design} \ {\rm of} \ {\rm Body}: {\rm Made} \ {\rm by} \ {\rm Aluminum} \ {\rm Die-Casting}$
- $\ast\,$ Powerful Motor for Stirring Quantities of up to 20Lit. (Water)
- * Ceramic Coated Plate-Chemical/Acid Resistance
- * Free Voltage (AC 80V ~ 240V, 50 / 60Hz) by SMPS
- * Included 1ea 3cm Magnetic Stirrer Bar

Specification

	"MS-20A"
Plate Size (mm)	180×180
Stirring Capacity (H2O)	Max. 20Lit.
Speed Range	Max. 1,500 rpm
Motor Rating	9/4W
Controller	Superior Linearity PWM Controller
Stirring Bar	3cm bar-included, Usable : up to 5cm bar
Permissible Ambient Temperature	Ambient +5 ∼ 50 °C
Permissible Relative Moisture	85%
Plate	Ceramic Coated Aluminum (Chemical / Acid Resistance)
Body	Powder Coated Aluminum
Dimension (w×d×h) (mm)	206×307×99
Weight (kg)	3.1 kg
Others	CE Certified PL(Product Liability) Insurance
Packing Size (w \times d \times h) (mm)	400×280×200
& Gross Weight	5 kg
Power Consumption	10 W
Power Supply & Cord/Plug*	Free Voltage (1 Phase, AC 80 ~ 240 V, 50 / 60 Hz)

* Other Specifications are available upon Customer's Request.

Ordering Information

DH.WMS03122 DH.WMS02122	Magnetic Stirrer, 180×180 mm, <code>230V, "MS-20A"</code> , with stir bar , <code>120V, "MS-20A"</code> , with stir bar
	Optional Accessories
DH.WHP503036	
	Φ12×L450mm, for " HP-20A/20D ", " MSH-20A/20D ", and " SMHS "
DH.WHP503012	(2) Holder / Clamp / Clip, "CL220", for Temp. Probe