


# InLab®605-10m

## dissolved oxygen



**Order number** 51340292

former order number -

### Specifications

measurement range 0 ... 200 % 0 ... 20 mg/L

compatible with SevenGo pro™ SG6 and  
SevenGo Duo pro™ SG68

temperature range 0 ... 60

temperature probe NTC 22 kΩ

shaft material PPS

shaft length 120 mm

shaft diameter 12 mm

connector

electrode head -

cable length 10.0 m

connector BNC/RCA (cinch)

### Storage

dry

### Special features

IP67

fixed cable

membrane is protected by a steel mesh

shaft made of glass fiber reinforced PPS

### Choice of accessories

Membrane kit: 51340293 (containing 3

membrane bodies and 25 mL of electrolyte)

Electrolyte: 51340294

Scale 1:1