


LE407**2-in-1****pH****Order number** 51340330

former order number -

Specifications

pH range 0 - 14

zero point (pH) 7.00 ± 0.25

temperature range 0 - 80 °C

temperature probe -

shaft material POM

membrane glass -

membrane resistance < 250 M Ω

reference system ARGENTHAL™

reference diaphragm Nylon

reference electrolyte 9811 paste

bridge electrolyte -

electrode head LE cap

cable length 1.0 M

connector BNC

cell constant (cm⁻¹) -

cell type -

indicating metal -

Storage

FRISCOLYT-B®

(51340053)

Special features

Unbreakable, chemically resistant POM shaft

Two nylon diaphragm

Scale 1:1