

Cat. No	Description	Articles
---------	-------------	----------

## DAIHAN-brand® Aluminum-case Stirrer Reaction Vessel Heating Mantles, "WHM", 450°C, with Certi. & Traceability For Reaction Vessel, Cylindrical-type, with Stirring & Temp. Controller, and Nickel Chrome Heating Element, 500 ~ 7,000ml

- with Certificate & Traceability : Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Full Choices of Vessel, Capacity, Housing and Control are Offered!
- **Nickel Chrome Heating Element - Coil Type** : Acid-Resistant, Long Life Time
- CE Certified
- **External Thermo Sensor is (K-type) available**, which Enables Precise Control of the Temperature of Liquid Sample in the Vessel when used together with Remote Digital Controller

Catalog No.		Flask Capacity (ml)	Flask Dia. $\phi$ (cm)	Watts (W)	Inside Depth (mm)	Mantle Outside $\phi \times H$ (mm)	Support Rod $\phi \times H$ (mm)
with 230V	with 120V						
DH.WHM722100	DH.WHM12210A	500	$\phi 11 \times H10$	240	72	$\phi 190 \times H165$	$\phi 13 \times H660$
DH.WHM722101	DH.WHM12210B	1,000	$\phi 11 \times H15$	350	107	$\phi 270 \times H200$	$\phi 13 \times H660$
DH.WHM722102	DH.WHM12210C	2,000	$\phi 13 \times H19$	450	142	$\phi 270 \times H235$	$\phi 13 \times H660$
DH.WHM722103	DH.WHM12210D	3,000	$\phi 13 \times H28$	550	177	$\phi 325 \times H270$	$\phi 13 \times H800$
DH.WHM722104	DH.WHM12210E	5,000	$\phi 16 \times H29$	800	177	$\phi 325 \times H270$	$\phi 13 \times H800$
DH.WHM722105	DH.WHM12210F	6,000	$\phi 21.5 \times H24$	800	177	$\phi 325 \times H270$	$\phi 13 \times H800$
DH.WHM722106	DH.WHM12210G	7,000	$\phi 21.5 \times H25$	850	177	$\phi 325 \times H270$	$\phi 13 \times H800$



(1) Analog Compact Controller &  
(2) Digital Feedback Controller  
(Optional, See P.90)



## DAIHAN-brand® Fabric-Housed Reaction Vessel Heating Mantles, "WHM", 0.5 ~ 7Lit, 450°C For Reaction Vessel, Cylindrical-type, with Nickel Chrome Heating Element, Option-Controller, with Certi. & Traceability

- with Certificate & Traceability : Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Ideal for Use with Magnetic Stirrer
- **Nickel Chrome Heating Element - Coil Type** : Acid-Resistant, Long Life Time
- Optional Remote Thermo Controllers are Available (Analog : WHM-C10A, Digital : WHM-C10D)
- **External Thermo Sensor is (K-type) available**, which Enables Precise Control of the Temperature of Liquid Sample in the Vessel when used together with Remote Digital Controller
- CE Certified

Catalog No.		Vessel Capacity (ml)	Vessel dia. $\phi \times$ height (cm)	Watts (W)	Inside Depth (mm)	Mantle Outside $\phi \times H$ (mm)
with 230V	with 120V					
DH.WHM12251	DH.WHM12261	500	$\phi 11 \times H10$	240	75	$\phi 160 \times H100$
DH.WHM12252	DH.WHM12262	1,000	$\phi 11 \times H15$	350	115	$\phi 160 \times H140$
DH.WHM12253	DH.WHM12263	2,000	$\phi 13 \times H19$	450	145	$\phi 185 \times H165$
DH.WHM12254	DH.WHM12264	3,000	$\phi 13 \times H28$	550	210	$\phi 185 \times H235$
DH.WHM12255	DH.WHM12265	4,000	$\phi 16 \times H26$	700	210	$\phi 220 \times H235$
DH.WHM12256	DH.WHM12266	5,000	$\phi 16 \times H29$	800	210	$\phi 220 \times H235$
DH.WHM12257	DH.WHM12267	6,000	$\phi 21.5 \times H24$	800	210	$\phi 290 \times H235$
DH.WHM12258	DH.WHM12268	7,000	$\phi 21.5 \times H25$	850	210	$\phi 290 \times H235$

● Option : Digital Feedback Controller or Analog Controller (Analog : WHM-C10A, Digital : WHM-C10D/ See P.90)