

DAIHAN-brand® Analog Aluminum-case Multi Heating Mantles, “WHM”, with 3 & 6 Places, 450°C, with Certi. & Traceability With Glassware Supporting Rod & Clamp, Ideal for Extraction · Reflux · Distillation Applications, Lower Profile, 250~1,000ml

- with Certificate & Traceability : Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Lower Profile for Space-saving Convenience
- Ideal for Repetitive Extraction (Kjeldahl, Soxhlet and other Extraction Procedure), Reflux and Distillation Laboratory Applications
- Analog Controller Built-in System
- Rods and Clamps are Offered as Standard
- Ideal for DURAN / Kimax / Pyrex & Iwaki-brand Round Flasks
- Highly-processed Aluminum Case Offers Robust Standing, Easy Handling
- **Nickel Chrome Heating Element - Coil Type : Acid-Resistant, Long Life Time**
- CE Certified
- with Glassware Supporting Rod & Clamp



Premium Analog Multi Heating Mantles

**Highly-processed Aluminum Case
Nickel Chrome Heating Element**



6 places Analog Multi Heating Mantle, 250ml, 300ml & 500ml, 450°C

Catalog No.		Flask Capacity (ml)	Flask		Watts(W)	Inside Depth (mm)	Dimension W×D×H (mm)
with 230V	with 120V		Max. Dia. φ (mm)	Bottom			
DH.WHM12291	DH.WHM12292	250/300	83 ~ 86	Flat	840	50	800×240×110
DH.WHM12293	DH.WHM12294	250/300	83 ~ 86	Round	840	50	800×240×110
DH.WHM12295	DH.WHM12296	500	102 ~ 105	Round	1440	65	950×270×130

3 places Analog Multi Heating Mantle, 500ml & 1,000ml, 450°C

Catalog No.		Flask Capacity (ml)	Flask		Watts(W)	Inside Depth (mm)	Dimension W×D×H (mm)
with 230V	with 120V		Max. Dia. φ (mm)	Bottom			
DH.WHM12391	DH.WHM12392	500	102 ~ 105	Round	720	65	550×270×130
DH.WHM12393	DH.WHM12394	1,000	130 ~ 133	Round	1050	75	690×340×170

※ Remarks : DU. : DURAN Ki. : Kimax Iw. : Iwaki Py. : Pyrex

Spare Support Rod & Clamp Set

- DH.WHM30500** Support Rod & Clamp Set, “MHR01”, for 3Places 500ml, 2×750-/1×500-mm Rod and 2×Clamp for Rod
- DH.WHM31000** , “MHR02”, for 3Places 1000ml, 2×750-/1×600-mm Rod and 2×Clamp for Rod
- DH.WHM60250** , “MHR03”, for 6Places 250-/300-ml, 2×750-/1×900-mm Rod and 2×Clamp for Rod
- DH.WHM60500** , “MHR04”, for 6Places 500ml, 2×750-/1×940-mm Rod and 2×Clamp for Rod