

Cat. No	Description	Articles
---------	-------------	----------

DAIHAN-brand® Heating Tapes, “WHM”, Standard-type, 250°C, 13 & 25mm Width

Catalog No.	13mm Width × Length(m)	Watts & Volt
DH.WHM12311	13mm × 0.6m	200W, 230V
DH.WHM12312	× 1.0m	200W, 230V
DH.WHM12313	× 1.5m	250W, 230V
DH.WHM12314	× 2.0m	300W, 230V
DH.WHM12315	× 2.5m	350W, 230V
DH.WHM12316	× 3.0m	400W, 230V
Catalog No.	25mm Width × Length(m)	Watts & Volt
DH.WHM12321	25mm × 0.6m	200W, 230V
DH.WHM12322	× 1.0m	200W, 230V
DH.WHM12323	× 1.5m	250W, 230V
DH.WHM12324	× 2.0m	300W, 230V
DH.WHM12325	× 2.5m	350W, 230V
DH.WHM12326	× 3.0m	400W, 230V
DH.WHM12327	× 4.0m	550W, 230V
DH.WHM12328	× 6.0m	700W, 230V

* For Glass Tubes Metal Pipet



Heating Tapes

DAIHAN-brand® Custom-Made Heating Mantles & Insulation

Bottom Valved Vessels or General-use, “WHM”, Up to 450°C ~ 650°C

We can Offer DAIHAN Custom Heating Mantles, for which you can Specify the Outside Diameter of the Vessel, Length of Mantle, Voltage, Number of Circuits and Length of Each Heating Zone, Operating Temperature, Unusual Ambient Conditions, Size & Location of Holes for View Ports, and Tubulations.

Please Submit a Sketch for Complicated Applications and Contact our International Sales Team (E-mail : daihan@daihansci.co.kr, Call Center : 080-008-3000) or our Distributor in your Country.

- Provides Uniform Heating and Insulation for Scrubbing and Fractionating Columns.
- Closes with a Zipper, with Eyelets and laces or Hooks and laces.
- Internal Operating Temperature is up to 450°C or 650°C

Examples for Customized Heating Mantles



Welcome! Custom-Made Heating Mantles!!