

SINCE 1979
FUYUE

Basic Series

Technical Specification:
 Material: Stainless Steel
 Class: F1、F2、M1
 Nominal Value: 1mg-50kg
 Density: 7.9g/cm³
 Susceptibility: 0.05-0.2
 Magnetization: F1 ≤ 25 μ T; F2 ≤ 80 μ T; M1 ≤ 250 μ T
 Package: Aluminum Box、Plastic Box
 Inlay: From UK. Approved by FDA and SGS



11-12

SINCE 1979
FUYUE

Patent Series

Technical Specification:
 Material: Stainless Steel
 Class: F1、F2、M1
 Nominal Value: 10g-20kg
 Density: 7.9-7.96g/cm³
 Susceptibility: ≤ 0.01
 Magnetization: F1 ≤ 25 μ T; F2 ≤ 80 μ T; M1 ≤ 250 μ T
 Package: Aluminum Box or Plastic Box
 Inlay: From UK. Approved by FDA and SGS



13-14



SINCE 1979
FUYUE