

SINCE 1979
FUYUE

High Grade Series

Technical Specification:

Material: High Grade Non-magnetic Stainless Steel

Class: E1, E2, F1

Nominal Value: 1mg-20kg

Density: $8.0 \pm 0.02 \text{ g/cm}^3$

Susceptibility: ≤ 0.005

Magnetization: $E1 \leq 2.5 \mu \text{ T}$; $E2 \leq 8 \mu \text{ T}$; $F1 \leq 25 \mu \text{ T}$

Package: Aluminum Box, Approved by FDA and SGS

Inlay: From UK. Approved by FDA and SGS



07-08

SINCE 1979
FUYUE

Standard Series

Technical Specification:

Material: Special Stainless Steel

Class: E2, F1

Nominal Value: 1mg-50kg

Density: $7.94-7.96 \text{ g/cm}^3$

Susceptibility: ≤ 0.01

Magnetization: $E2 \leq 8 \mu \text{ T}$; $F1 \leq 25 \mu \text{ T}$

Package: Aluminum Box, Approved by FDA and SGS

Inlay: From UK. Approved by FDA and SGS



10kg-20kg



09-10