

SINCE 1979  
**FUYUE**

## Rectangular Series

Technical Specification:  
Material: Stainless Steel  
Class: F1, F2, M1  
Nominal Value: 5kg-20kg  
Density: 7.92g/cm<sup>3</sup>  
Susceptibility: <0.1  
Package: Aluminum Box, Paper Box  
Inlay: From UK. Approved by FDA and SGS




5kg



10kg



20kg



15-16

SINCE 1979  
**FUYUE**

## Chromed Steel Series

Technical Specification:  
Material: Iron/Steel  
Class: M1  
Nominal Value: 1g-20kg  
Density: 7.85 g/cm<sup>3</sup>  
Magnetization: M1 ≤ 250 μ T  
Surface: Chromed  
Package: Plastic Box  
Inlay: From UK. Approved by FDA and SGS



17-18