


InLab[®] Expert Pt1000

pH Ref °C ISM



For general applications in the lab and outdoor;
robust and maintenance-free

Scale 1:1



Sensor type combined pH electrode

Order number 51343105

Specifications

pH range 0...14

zero point (pH) 7.00 ± 0.25

temperature range 0...100 °C

temperature probe Pt1000

shaft material PEEK

membrane glass U

membrane resistance $< 250 \text{ M}\Omega$

reference system ARGENTHAL[™]

reference diaphragm double open junction

reference electrolyte XEROLYT[®] polymer

bridge electrolyte -

electrode head MultiPin[™]

cable length -

connector -

cell constant (cm⁻¹) -

cell type -

indicating metal -

Storage

InLab storage solution (30111142)

or FRISCOLYT-B[®] (51340053)