




Scale 1:1



Sensor type combined pH electrode

Order number 59902917

Specifications

pH range 1...11

zero point (pH) 7.00 ± 0.25

temperature range 0...80 °C

temperature probe -

shaft material glass

membrane glass U

membrane resistance < 250 MΩ

reference system SteadyForce™

reference diaphragm double ceramic

reference electrolyte DPA gel

bridge electrolyte -

electrode head S7

cable length -

connector -

cell constant (cm⁻¹) -

cell type -

indicating metal -

Storage

InLab storage solution (30111142)
or 3 mol/L KCl (51350072)