


InLab[®] Hydrofluoric

For hydrofluoric acid containing samples

pH Ref °C ISM



Scale 1:1



Sensor type pH half-cell

Order number 51343176

Specifications

pH range 1...11

zero point (pH) 7.00 ± 0.25

temperature range 0...100 °C

temperature probe -

shaft material glass

membrane glass HF

membrane resistance < 100 MΩ

reference system ARGENTHAL™

reference diaphragm ceramic

reference electrolyte 3 mol/L KCl

bridge electrolyte -

electrode head S7

cable length -

connector -

cell constant (cm⁻¹) -

cell type -

indicating metal -

Storage

InLab storage solution (30111142)

or 3 mol/L KCl (51350072)