


InLab[®] Pure

pH Ref °C ISM



For ion-weak-samples like pure water, deionized water

Scale 1:1



Sensor type combined pH electrode

Order number 30248112

Specifications

pH range 1...11

zero point (pH) 7.00 ± 0.25

temperature range 0...80 °C

temperature probe -

shaft material glass

membrane glass Lot

membrane resistance < 50 MΩ

reference system ARGENTHAL™

reference diaphragm fixed glass sleeve

reference electrolyte 3 mol/L KCl

bridge electrolyte -

electrode head S7

cable length -

connector -

cell constant (cm⁻¹) -

cell type -

indicating metal -

Storage

InLab storage solution (30111142)

or 3 mol/L KCl (51350072)