



Scale 1:1

**Order number** 30014092  
 former order number -

**Specifications**  
 pH range -  
 zero point (pH) -

temperature range 0 - 100 °C  
 temperature probe NTC 30 kΩ

shaft material epoxy  
 membrane glass -  
 membrane resistance -

reference system -  
 reference diaphragm -  
 reference electrolyte -

bridge electrolyte -

electrode head -  
 cable length 1.2 m  
 connector Mini-DIN

cell constant (cm<sup>-1</sup>) 0.57  
 cell type 4 graphite poles

indicating metal -

**Storage** dry

**Additional specifications**  
 measuring range 0.01 - 1000 mS/cm  
 compatible with appropriate METTLER TOLEDO meters with Intelligent Sensor Management (ISM) functionality