



Scale 1:1

Order number 51302256

former order number -

**Specifications**

pH range -

zero point (pH) -

temperature range 0 - 100 °C

temperature probe NTC 30 kΩ

shaft material glass

membrane glass -

membrane resistance -

reference system -

reference diaphragm -

reference electrolyte -

bridge electrolyte -

electrode head -

cable length 1.0 m

connector Mini-DIN

cell constant (cm<sup>-1</sup>) 0.80

cell type 4 platinum poles

indicating metal -

**Storage** dry**Additional specifications**

measuring range 0.01 - 500 mS/cm