



Scale 1:1

Order number 30014160

former order number -

Specifications

pH range -

zero point (pH) -

temperature range 0 - 100 °C

temperature probe NTC 30 kΩ

shaft material glass

membrane glass -

membrane resistance -

reference system -

reference diaphragm -

reference electrolyte -

bridge electrolyte -

electrode head -

cable length 1.0 m

connector Mini-DIN

cell constant (cm⁻¹) 0.1

cell type 2 platinum poles

indicating metal -

Storage dry**Additional specifications**

measuring range 0.1 - 500 μS/cm

especially for the measurement of conductivity
in bio-ethanolcompatible with METTLER TOLEDO SevenCom-
pact™ conductivity meter S230, SevenMulti™,
FiveEasy™ and FiveGo™.

platinum 2-plate conductivity probe

integrated temperature sensor

fixed cable