



Order number 51344128
 former order number

Specifications
 pH range -
 zero point (pH) -

temperature range 0 - 100 °C
 temperature probe NTC 30 kΩ

shaft material stainless steel
 membrane glass -
 membrane resistance -

reference system -
 reference diaphragm -
 reference electrolyte -

bridge electrolyte -

electrode head -
 cable length 5.0 m
 connector Mini-DIN

cell constant (cm⁻¹) 0.105
 cell type 2 steel poles

indicating metal -

Storage dry

Additional specifications
 measuring range 0.001 - 500 μS/cm
 IP67

Scale 1:1