



Scale 1:1

Order number 51344118

former order number

Specifications

pH range -
 zero point (pH) -

temperature range 0 - 100 °C

temperature probe NTC 30 kΩ

shaft material stainless steel

membrane glass -

membrane resistance -

reference system -

reference diaphragm -

reference electrolyte -

bridge electrolyte -

electrode head -

cable length 5.0 m

connector Mini-DIN (IP67)

cell constant (cm⁻¹) 0.105

cell type 2 steel poles

indicating metal -

Storage dry

Additional specifications

measuring range 0.001 - 500 µS/cm

compatible with appropriate METTLER TOLEDO meters with Intelligent Sensor Management (ISM) functionality

IP67