

## Accessories Series



stainless steel sling



stainless steel box



steel painting box



heavy capacity aluminum box



Brush



Dust Remover



Non-Magnetic Tweezer



Handle



Glove



Cloth



## Accessories Series

## Weight Sequence

Nominal Value	1mg-500mg	1mg-100g	1mg-200g	1mg-500g	1mg-1kg	1mg-2kg	1mg-5kg	1mg-10kg	1kg-5kg
1mg	1	1	1	1	1	1	1	1	---
2mg	2	2	2	2	2	2	2	2	---
5mg	1	1	1	1	1	1	1	1	---
10mg	1	1	1	1	1	1	1	1	---
20mg	2	2	2	2	2	2	2	2	---
50mg	1	1	1	1	1	1	1	1	---
100mg	1	1	1	1	1	1	1	1	---
200mg	2	2	2	2	2	2	2	2	---
500mg	1	1	1	1	1	1	1	1	---
1g	---	1	1	1	1	1	1	1	---
2g	---	2	2	2	2	2	2	2	---
5g	---	1	1	1	1	1	1	1	---
10g	---	1	1	1	1	1	1	1	---
20g	---	2	2	2	2	2	2	2	---
50g	---	1	1	1	1	1	1	1	---
100g	---	1	1	1	1	1	1	1	---
200g	---	---	2	2	2	2	2	2	---
500g	---	---	---	1	1	1	1	1	---
1kg	---	---	---	---	1	1	1	1	1
2kg	---	---	---	---	---	2	2	2	2
5kg	---	---	---	---	---	---	1	1	1
10kg	---	---	---	---	---	---	---	1	---
20kg	---	---	---	---	---	---	---	---	---
Total	12	21	23	24	25	27	28	29	4

## OIML MPE (±mg)

Weight class	E1	E2	F1	F2	M1	M2
2000kg			10000	30000	100000	300000
1000kg		1600	5000	16000	50000	160000
500kg		800	2500	8000	25000	80000
200kg		300	1000	3000	10000	30000
100kg		160	500	1600	5000	16000
50kg	25	80	250	800	2500	8000
20kg	10	30	100	300	1000	3000
10kg	5	16	50	160	500	1600
5kg	2.5	8	25	80	250	800
2kg	1	3	10	30	100	300
1kg	0.5	1.6	5	16	50	160
500g	0.25	0.8	2.5	8	25	80
200g	0.1	0.3	1	3	10	30
100g	0.05	0.16	0.5	1.6	5	16
50g	0.03	0.1	0.3	1	3	10
20g	0.025	0.08	0.25	0.8	2.5	8
10g	0.02	0.06	0.2	0.6	2	6
5g	0.016	0.05	0.16	0.5	1.6	5
2g	0.012	0.04	0.12	0.4	1.2	4
1g	0.01	0.03	0.1	0.3	1	3
500mg	0.008	0.025	0.08	0.25	0.8	2.5
200mg	0.006	0.020	0.06	0.2	0.6	2
100mg	0.005	0.016	0.05	0.16	0.5	1.6
50mg	0.004	0.012	0.04	0.12	0.4	1.2
20mg	0.003	0.010	0.03	0.1	0.3	1.0
10mg	0.003	0.008	0.025	0.08	0.25	0.8
5mg	0.003	0.006	0.02	0.06	0.2	0.6
2mg	0.003	0.006	0.02	0.06	0.2	0.6
1mg	0.003	0.006	0.02	0.06	0.2	0.6

## ASTM MPE (±mg)

Weight class	Class 0	Class 1	Class 2	Class 3	Class 4	Class 5
50kg	63	125	250	500	1000	2500
20kg	25	50	100	200	400	1000
10kg	13	25	50	100	200	500
5kg	6	12	25	50	100	250
2kg	2.5	5	10	20	40	100
1kg	1.3	2.5	5	10	20	50
500g	0.6	1.2	2.5	5	10	30
200g	0.25	0.5	1	2	4	15
100g	0.13	0.25	0.5	1	2	9
50g	0.06	0.12	0.25	0.6	1.2	5.6
20g	0.037	0.074	0.1	0.35	0.7	3
10g	0.025	0.05	0.074	0.25	0.5	2
5g	0.017	0.034	0.054	0.18	0.36	1.3
2g	0.017	0.034	0.054	0.13	0.26	0.75
1g	0.017	0.034	0.054	0.1	0.2	0.5
500mg	0.005	0.01	0.025	0.08	0.16	0.38
200mg	0.005	0.01	0.025	0.06	0.12	0.26
100mg	0.005	0.01	0.025	0.05	0.1	0.2
50mg	0.005	0.01	0.014	0.042	0.085	0.16
20mg	0.005	0.01	0.014	0.035	0.07	0.12
10mg	0.005	0.01	0.014	0.03	0.06	0.1
5mg	0.005	0.01	0.014	0.028	0.055	0.08
2mg	0.005	0.01	0.014	0.025	0.05	0.06
1mg	0.005	0.01	0.014	0.025	0.05	0.05