## B-500BPh





Plan-achromatic objectives with mechanical tube lenght 160 mm, made by following objectives: -) Plan-achromatic 10XPh, N.A. 0.25, W.D. 1,30 mm -) Plan-achromatic 20XPh, N.A. 0.40, W.D. 0,55 mm -) Plan-achromatic 40XPh, N.A. 0.65, W.D. 0,25 mm -) Plan-achromatic 100XPh, N.A. 1,25, W.D. 0,10 mm (oil immersion) All objectives are treated with an anti-fungus treatment. Condenser Phase contrast condenser (for 10x, 20x, 40x, 100x objectives) and brightfield. Condenser for brightfield, swing-out N.A. 1,2. Dimensions HEIGHT: 420 mm WIDTH: 250 mm DEPTH: 290 mm WEIGHT: 8 kg

## Accessories Centering Telescope included. Instruction manual and dust cover included.